

3: Perforated Tray for 2.5 gal and SS 2.5/3.5 gal Cleaners **Cole-Parmer Cat No: EW-08853-42**

Features:

- Made of 304 stainless steel, these insert trays can be used to support items that require indirect cleaning.
- Use perforated or mesh trays to support items in the tank for direct cleaning.
- Never use ultrasonic cleaners with flammable solvents. A cleaning solution must be used with these cleaners.
- Use an insert tray or beaker positioning cover and beaker with these ultrasonic cleaners.
- Ultrasonic activity is greatest starting about an inch from the bottom of the tank.



Specifications:

Description	Perforated Tray for 2.5 gal and SS 2.5/3.5 gal Cleaners
Origin	U.S.A